

T E C H N O P A T H


MULTICHEM IA: PROVIDING TEST CONSOLIDATION FOR IMMUNOASSAY QC IN A LIQUID STABLE FORMAT

Using fresh frozen, human based plasma and a proprietary formulation, TECHNOPATH introduce a liquid stable, tri-level, immunoassay control, with unsurpassed stability, performance and analyte consolidation. Including Thyroids, Fertility, TDM's, Haematinics, Tumor Markers and Cardiac Markers in the one control, Multichem IA maximises the cost savings resulting from reduced labour and inventory, coupled with increased efficiency on your automated platform. Its unique formulation means that Multichem IA behaves just like your patient samples and thereby provides you with added confidence you require in your results and eliminates the inefficiencies and costs associated with reruns due to poor QC performance. MULTICHEM IA provides the essential features to fulfil the needs of today's laboratory, where platform consolidation and premium performance are vital components to providing optimal patient care and cost efficiency.

FEATURES:

- Proprietary Liquid Stable Formulation Providing Unsurpassed Stability and Performance.
- Tri-level control 12 x 5ml Kit Size.
- Dropper System for Convenient and Optimised Usage.
- 10 day Open vial stability at 2-8C.
- 30 Month Closed Vial Stability at -20C.
- Human based material providing excellent accuracy and precision.
- Kits are barrier wrapped for added protection and stability during shipment and storage.
- Customised and/or OEM lot numbers available depending on volume and analyte range.

TECHNOPATH, Fort Henry Business Park, Ballina, Co. Tipperary.
Ph: +353-61-335844 U.K.: +44 (0)28 30833808 Fax: +353-61-203034
Email: info@techno-path.com www.techno-path.com



T E C H N O P A T H

MULTICHEM IA: IMMUNOASSAY ANALYTE LIST

Cancer Markers	Adrenal Pituitary	Renal
Alpha Fetoprotein	Cortisol	Aldosterone*
Carcinogenic Embryonic Anti-gen	Adrenocorticotrophic Hormone (ACTH)*	Renin*
Prostate Specific Antigen, Total		
Prostate Specific Antigen, Free		
CA 125		
CA 19-9		
CA 15-3		
Cardiac	Drugs	Bone Metabolism
CK-MB	Amikacin	Parathyroid H (Intact PTH)
Troponin I	Carbamazepine	Human Growth Hormone
Troponin T	Caffeine	Insulin-Like Growth factor (IGF1)
Pro-BNP/BNP*	Chloramphenicol	
UltraSensitive CRP	Digoxin	
Myoglobin	Disopyramide	
Allergy	Ethosuximide	Reproductive/ Fertility
IgE*	Gentamicin	Estradiol
	Ibuprofen	Human Chorionic Gonadotropin
	Lidocaine	Progesterone
	Lithium	Prolactin
	N-Acetyl procainamide	Follicle Stimulating Hormone
	Paracetamol	Luteinising Hormone
	Phenobarbital	Testosterone
	Procainamide	DHEA Sulfate
	Primidone	Sex Hormone Binding Globulin
	Phenytoin	17 Hydroxyprogesterone
	Quinidine	
	Salicylate	Diabetes
	Theophylline	Insulin
	Tobramycin	C-Peptide
	Valproic Acid	
	Vancomycin	Anaemia
		Ferritin
		Folate
		Vitamin B12

Note* Available in specific lots only.