

T E C H N O P A T H


MULTICHEM U: PROVIDING TEST CONSOLIDATION FOR URINARY CHEMISTRY QC IN A LIQUID STABLE FORMAT

Using a proprietary formulation, TECHNOPATH introduce a liquid stable, human based, bi-level, urine chemistry control with unsurpassed stability, performance and analyte consolidation. MULTICHEM U provides the added confidence you require in your results and eliminates the inefficiencies and costs associated with reruns due to poor QC performance. MULTICHEM U provides the essential features to fulfil the needs of today's laboratory, where platform consolidation and premium performance are vital components to providing optimal patient care and cost efficiency.

FEATURES:

- Proprietary liquid stable formulation providing unsurpassed stability and performance
- 2 Levels; Available in 15 x 10 ml Kit Size
- Optional dropper system for convenient and optimised usage
- 30 day open vial stability at 2-8C
- 2 year closed vial stability at 2-8C
- Human based analytes including Microalbumin provides excellent accuracy and reproducibility
- Kits are barrier wrapped for added protection and stability during shipment and storage.
- Customised and/or OEM lot numbers available depending on volume and analyte range.

TECHNOPATH, Unit 3, Fort Henry Business Park, Ballina, Co. Tipperary.
Ph: +353-61-335844 U.K.: +44 (0)28 30833808 Fax: +353-61-203034
Email: info@techno-path.com www.techno-path.com



T E C H N O P A T H

MULTICHEM-U URINE CHEMISTRY ANALYTE LIST

Amylase

Calcium

Chloride

Cortisol

Creatinine

Glucose

Magnesium

Microalbumin

Osmolality

Phosphate

Potassium

Sodium

Urea

Uric acid

Urinary Protein

HCG *

ALA *

VMA *

HVA *

5-HIAA *

Note: * Available in specific lots only.

TECHNOPATH, Unit 3, Fort Henry Business Park, Ballina, Co. Tipperary.
Ph: +353-61-335844 U.K.: +44 (0)28 30833808 Fax: +353-61-203034
Email: info@techno-path.com www.techno-path.com